




STERILIZATION
SYSTEMS

Mini Bulk and Bulk Size
STEAM STERILIZERS

LARGEST SELECTION OF SIZES AVAILABLE ON THE MARKET




STERILIZATION
SYSTEMS

125 James Way, Southampton, PA. USA 18966
Ph. 215.355.9100 | F. 215.357.4000
etcSterilization.com | © 2020 ETC.




ETC GIVES YOU MORE CHOICES


ETC offers more sterilizer chamber sizes than the competition, multiple door designs that fit any space, and sectional autoclaves that are ideal for hard to reach areas.



A



B



C

THREE UNIQUE DOOR DESIGNS

- A. Vertical Rising Door Design
- B. Hinged Door Design
- C. Horizontal Sliding Door Design

LARGEST SELECTION OF SIZES

ETC offers over 60 cross sectional sizes ranging from mini bulk (36" W x 42" H) to bulk sized units (48" W x 84" H, 72" W x 84" H) plus

CHOICE OF CONTROL SYSTEMS

- IDEC Microsmart PLC
- Allen Bradley Compact Logix

EASY GLIDE LOADING CART

This heavy duty loading cart eliminates the need for the typical loading cart/loading carriage/trolley combination for floor mounted autoclaves



SECTIONAL AUTOCLAVES

ETC's process for installing Bulk Size autoclaves in basements, on upper floors and in other challenging facility locations is based on the practice of installing the sterilizer in modular segments or sections.

This unique procedure eliminates the need for rigging full size units through openings in facility walls or ceilings, costly facility modifications or any disruption to everyday business activities surrounding the installation. Each section is delivered to the final destination with no additional modifications.

